



Telephone: (650) 697-3600

GOAT ANTI-14-3-3 TAU / YWHAQ ANTIBODY

SKU: EB08137



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

tyrosine 3/tryptophan 5 -monooxygenase activation protein, theta polypeptide|protein tau|14-3-3 protein theta|14-3-3 protein T-cell|HS1|1C5|14-3-3|tyrosine 3-monooxygenase/tryptophan 5-monooxygenase

activation protein, theta polypeptide|14-3-3 protein tau|YWHAQ Names

Accession

Alias

ID

NP_006817.1

Blocking

EBP08137 **Peptide**

Peptide with sequence C-DDRKQTIDNSQ, from the internal region of the protein sequence according to **Immunogen**

NP 006817.1.

Peptide

C-DDRKQTIDNSQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog **Species**

Reactive

Human **Species**

Human

10971 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Antibody detection limit dilution 1:4000.

Limit

Western Approx 28kDa band observed in Human Brain (Hippocampus) lysates (calculated MW of 27.8kDa according to

Blot NP_006817.1). Recommended concentration: 0.03-0.1µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa