

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-A4GNT / ALPHA4GNT ANTIBODY**

**SKU:** EB09734



Telephone: (650) 697-3600



250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / **Alias** 

alpha4GnT|MGC149493|alpha-1, 4-N-acetylglucosaminyltransferase|A4GNT

Names

ID

**Accession** 

NP\_057245.1

**Blocking** 

EBP09734 **Peptide** 

Peptide with sequence C-PRTYRDLIKGPEGS, from the C Terminus of the protein sequence according to **Immunogen** 

NP 057245.1.

**Peptide** 

C-PRTYRDLIKGPEGS Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Dog, Cow **Species** 

Reactive

Human **Species** 

Human

51146 **Gene ID** 

**Product** 

Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:8000.

Limit

Western Approx 40kDa band observed in lysates of cell line K562 (calculated MW of 39.5kDa according to

**Blot** NP\_057245.1). Recommended concentration: 1-3µg/ml.

**Application** 

Pep-ELISA, WB Type

## **GALLERY IMAGES**



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa