



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-AADAT ANTIBODY

SKU: EB08259



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

kynurenine aminotransferase II|alpha-aminoadipate aminotransferase|L-kynurenine/alpha-aminoadipate

Alias aminotransferase|L kynurenine/alpha aminoadipate aminotransferase|KATII|KAT2|aminoadipate

Names aminotransferase|AADAT

Accession

NP_057312.1; NP_872603.1

Blocking

EBP08259 **Peptide**

Peptide with sequence C-KPEDAKNPQKNTPK, from the internal region of the protein sequence according to **Immunogen**

NP_057312.1; NP_872603.1.

Product

Both reported variants represent identical protein (NP 057312.1and NP 872603.1). **Comments**

Peptide

Sequence

C-KPEDAKNPQKNTPK

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Reactive **Species**

Human

51166

Human

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western Approx 45kDa band observed in Human Liver lysates (calculated MW of 47.4kDa according to NP_057312.1).

Recommended concentration: 1-3µg/ml. Blot

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa