

GOAT ANTI-ABCA4 (AA1336-49) ANTIBODY

SKU: EB09205



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	retina-specific ABC transporter photoreceptor rim protein OTTHUMP00000012366 ATP-binding transporter, retina-specific ATP-binding cassette, sub-family A member 4 ATP binding cassette transporter STGD1 STGD RP19 RMP FLJ17534 FFM DKFZp781N1972 CORD3 ARMD2 ABCR ABC10 ATP-binding cassette, sub-family A (ABC1), member 4 ABCA4
Accession ID	NP_000341.2
Blocking Peptide	EBP09205
Immunogen	Peptide with sequence EPECPGPQLNTGTQ, from the internal region of the protein sequence according to NP_000341.2.
Peptide Sequence	EPECPGPQLNTGTQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	24
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)