

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ABCA4 (AA1336-49) ANTIBODY

SKU: EB09205

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

retina-specific ABC transporter|photoreceptor rim protein|OTTHUMP00000012366|ATP-binding transporter, retina-

specific|ATP-binding cassette, sub-family A member 4|ATP binding cassette

Alias transporter|STGD1|STGD|RP19|RMP|FLJ17534|FFM|DKFZp781N1972|CORD3|ARMD2|ABCR|ABC10|ATP-binding **Names**

cassette, sub-family A (ABC1), member 4|ABCA4

Accession

NP 000341.2 ID

Blocking

EBP09205 **Peptide**

Peptide with sequence EPECPGPQLNTGTQ, from the internal region of the protein sequence according to Immunogen

NP 000341.2.

Peptide

EPECPGPQLNTGTQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

24

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:2000.

Limit

Western Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this Blot product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.

Application

Type

Pep-ELISA