

GOAT ANTI-ABCB8 ANTIBODY

SKU: EB08731



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	mitochondrial ABC protein ATP-binding cassette, sub-family B, member 8 MABC1 M-ABC1 EST328128 ATP-binding cassette, sub-family B (MDR/TAP), member 8 ABCB8
Accession ID	NP_009119.2
Blocking Peptide	EBP08731
Immunogen	Peptide with sequence C-KKPEGPRSHQHK, from the C Terminus of the protein sequence according to NP_009119.2.
Peptide Sequence	C-KKPEGPRSHQHK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	11194
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)