

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ABCB8 ANTIBODY

SKU: EB08731

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

mitochondrial ABC protein|ATP-binding cassette, sub-family B, member 8|MABC1|M-ABC1|EST328128|ATP-binding **Alias** cassette, sub-family B (MDR/TAP), member 8|ABCB8

Names

Accession NP 009119.2 ID

Blocking EBP08731

Peptide

Peptide with sequence C-KKPEGPRSHQHK, from the C Terminus of the protein sequence according to Immunogen NP 009119.2.

Peptide

C-KKPEGPRSHQHK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Refrigerated Instructions

Predicted Human **Species**

Human 11194

Grade

Gene ID Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection Antibody detection limit dilution 1:128000. Limit

Application Pep-ELISA **Type**