



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ABCC1 ANTIBODY

SKU: EB08734







250kDa 150kDa 100kDa 75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / **Alias**

multiple drug resistance-associated protein|multiple drug resistance protein 1|multidrug resistance protein|leukotriene C(4) transporter|LTC4 transporter|ATP-binding cassette, sub-family C, member 1|MRP1|MRP|GS-X|DKFZp781G125|DKFZp686N04233|ABCC|ABC29|ATP-binding cassette, sub-family C

Names

(CFTR/MRP), member 1|ABCC1

Accession

NP_004987.2; NP_063915.2; NP_063953.2; NP_063954.2; NP_063955.2

Blocking

EBP08734 **Peptide**

Peptide with sequence C-HQSDLKVDENQKAYY, from the internal region of the protein sequence according to **Immunogen**

NP_004987.2; NP_063915.2; NP_063953.2; NP_063954.2; NP_063955.2.

Product This antibody is expected to recognize all reported isoforms (NP_004987.2; NP_063915.2; NP_063953.2;

NP_063954.2; NP_063955.2). Comments

Peptide

Sequence

C-HQSDLKVDENQKAYY

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog, Pig, Cow **Species**

Reactive

Species

Human Gene ID

4363

Human

Mouse Gene ID

17250

Rat Gene ID 24565

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:128000.

Limit

Western Blot

Approx 150kDa band observed in lysates of cell line A549 (calculated MW of 159kDa according to NP_063954.2). Recommended concentration:0.05-0.2µg/ml.

Application

Type

Pep-ELISA, WB



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa