


GOAT ANTI-ABCD2 (AA460-473) ANTIBODY

SKU: EB11507



250kDa
150kDa
100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	hALDR ATP-binding cassette, sub-family D (ALD), member 2 ATP-binding cassette sub-family D member 2 ALDRP ALDR ALDL1 adrenoleukodystrophy-related protein adrenoleukodystrophy-like 1 ABCD2 ABC39
Accession ID	NP_005155.1
Blocking Peptide	EBP11507
Immunogen	Peptide with sequence C-PLSDTLAIGKVID, from the internal region of the protein sequence according to NP_005155.1.
Product Comments	No cross-reactivity expected with ABCD1.
Peptide Sequence	C-PLSDTLAIGKVID
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	225
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 80kDa band observed in Human Liver lysates (calculated MW of 83.2kDa according to NP_005155.1). Recommended concentration: 0.5-2µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

