



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ABCF3 ANTIBODY

SKU: EB08739

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias FLJ11198|EST201864|ATP-binding cassette, sub-family F (GCN20), member 3|ABCF3

Names

Accession NP 060828.2 ID

Blocking

EBP08739 **Peptide**

Immunogen NP_060828.2. Peptide with sequence C-KESRLESSGKNKSYD, from the internal region of the protein sequence according to

Peptide C-KESRLESSGKNKSYD Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Refrigerated

Instructions

Predicted

Human, Mouse, Rat, Dog, Cow **Species**

Human 55324

Gene ID

Mouse 27406 **Gene ID**

Rat Gene ID 287982

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Application Pep-ELISA

Type