

GOAT ANTI-ABCF3 ANTIBODY

SKU: EB08739



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FLJ11198 EST201864 ATP-binding cassette, sub-family F (GCN20), member 3 ABCF3
Accession ID	NP_060828.2
Blocking Peptide	EBP08739
Immunogen	Peptide with sequence C-KESRLESSGKNKSYD, from the internal region of the protein sequence according to NP_060828.2.
Peptide Sequence	C-KESRLESSGKNKSYD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	55324
Mouse Gene ID	27406
Rat Gene ID	287982
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)