

# GOAT ANTI-ABR ANTIBODY

**SKU:** EB06101



## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                                                                                |
|------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                                                                                         |
| <b>Synonym /</b>             |                                                                                                                                                                                                                                                                                                |
| <b>Alias</b>                 | ABR active BCR-related gene MDB FLJ45954 active breakpoint cluster region-related protein                                                                                                                                                                                                      |
| <b>Names</b>                 |                                                                                                                                                                                                                                                                                                |
| <b>Accession ID</b>          | NP_068781.2; NP_001083.2; NP_001153218.1                                                                                                                                                                                                                                                       |
| <b>Immunogen</b>             | Peptide with sequence C-ELKRNTLYFSTDV, from the C Terminus of the protein sequence according to NP_068781.2; NP_001083.2; NP_001153218.1.                                                                                                                                                      |
| <b>Product</b>               | This antibody is expected to recognise all reported human isoforms of this protein (as represented by                                                                                                                                                                                          |
| <b>Comments</b>              | NP_068781.2; NP_001083.2 and NP_001153218.1).                                                                                                                                                                                                                                                  |
| <b>Peptide Sequence</b>      | C-ELKRNTLYFSTDV                                                                                                                                                                                                                                                                                |
| <b>Shipping Instructions</b> | Refrigerated                                                                                                                                                                                                                                                                                   |
| <b>Predicted Species</b>     | Human, Mouse, Rat, Dog, Cow                                                                                                                                                                                                                                                                    |
| <b>Human Gene ID</b>         | 29                                                                                                                                                                                                                                                                                             |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>                                                        |
| <b>ELISA</b>                 |                                                                                                                                                                                                                                                                                                |
| <b>Detection Limit</b>       | Antibody detection limit dilution 1:32000.                                                                                                                                                                                                                                                     |
| <b>Western Blot</b>          | Preliminary experiments in Human Brain lysates and in various cell line lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? |
| <b>Application Type</b>      | Pep-ELISA                                                                                                                                                                                                                                                                                      |