Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-ACACB (AA146-159) ANTIBODY**

**SKU:** EB11134

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

ACACB| HACC275| acetyl-Coenzyme A carboxylase beta| acetyl-CoA carboxylase beta| acetyl-CoA carboxylase 2| **Alias** 

ACC-beta| ACCB| ACC2

**Names** 

Accession NP 001084.3

**Blocking** 

EBP11134 **Peptide** 

Peptide with sequence C-SPSKEDKKQANIKR, from the internal reion (near N Terminus) of the protein sequence Immunogen

according to NP 001084.3.

**Peptide** 

C-SPSKEDKKQANIKR Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human **Species** 

Human

32 Gene ID

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:128000.

Limit

Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the Western

calculated size of 277kDa according to NP\_001084.3. Therefore we cannot recommend an optimal concentration Blot

and the antibody is an aspiring product.

Application

Pep-ELISA **Type**