

GOAT ANTI-ACACB (AA146-159) ANTIBODY

SKU: EB11134



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	ACACB HACC275 acetyl-Coenzyme A carboxylase beta acetyl-CoA carboxylase beta acetyl-CoA carboxylase 2 ACC-beta ACCB ACC2
Accession ID	NP_001084.3
Blocking Peptide	EBP11134
Immunogen	Peptide with sequence C-SPSKEDKKQANIKR, from the internal reion (near N Terminus) of the protein sequence according to NP_001084.3.
Peptide Sequence	C-SPSKEDKKQANIKR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	32
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 277kDa according to NP_001084.3. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)