

# GOAT ANTI-ACACB ANTIBODY

SKU: EB07748



## SPECIFICATIONS

<b>Formulation</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Synonym / Alias Names</b>	acetyl-CoA carboxylase 2 HACC275 ACCB acetyl-Coenzyme A carboxylase beta ACACB
<b>Accession ID</b>	NP_001084.2
<b>Blocking Peptide</b>	EBP07748
<b>Immunogen</b>	Peptide with sequence C-QLGEPDLSKDRKD, from the internal region of the protein sequence according to NP_001084.2.
<b>Peptide Sequence</b>	C-QLGEPDLSKDRKD
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Human Gene ID</b>	32
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:32000.
<b>Western Blot</b>	Preliminary experiments gave an approx 28kDa band in Human Heart lysates after 0.5µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 277kDa according to NP_001084.2. The 28kDa band was successfully blocked by incubation with the immunizing peptide.
<b>Application Type</b>	Pep-ELISA

## DOCUMENTS

- [Data Sheet](#)