

GOAT ANTI-ACE2 (INTERNAL) ANTIBODY

SKU: EB13086



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Angiotensin-converting enzyme homolog antibody Angiotensin I converting enzyme 2 antibody Angiotensin I Converting Enzyme (peptidyl dipeptidase A) 2 antibody Angiotensin converting enzyme like protein antibody Angiotensin converting enzyme homolog antibody Angiotensin converting enzyme 2 antibody ACEH antibody ACE2_HUMAN antibody ACE2 antibody ACE-related carboxypeptidase antibody ACE related carboxypeptidase antibody ACE 2 antibody
Accession ID	NP_001358344.1
Blocking Peptide	EBP13086
Immunogen	Peptide with sequence PROPRIETARY-ACE2, from the internal region of the protein sequence according to NP_001358344.1.
Peptide Sequence	PROPRIETARY-ACE2
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	59272
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 100kDa band was also observed on full-length recombinant Human ACE2 protein (calculated MW of 92.4kDa according to NP_001358344.1). Recommended concentration 1µg/ml. Primary incubation 1 hour at room temperature.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

