


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ACE2 (INTERNAL) ANTIBODY

SKU: EB13086







250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Angiotensin-converting enzyme homolog antibody|Angiotensin I converting enzyme 2 antibody|Angiotensin I

Synonym /

Converting Enzyme (peptidyl dipeptidase A) 2 antibody|Angiotensin converting enzyme like protein

antibody|Angiotensin converting enzyme homolog antibody|Angiotensin converting enzyme 2 antibody|ACEH Alias **Names** antibody|ACE2_HUMAN antibody|ACE2 antibody|ACE-related carboxypeptidase antibody|ACE related

carboxypeptidase antibody|ACE 2 antibody

Accession

NP_001358344.1

Blocking

EBP13086 **Peptide**

Peptide with sequence PROPRIETARY-ACE2, from the internal region of the protein sequence according to **Immunogen**

NP_001358344.1.

Peptide

Sequence

PROPRIETARY-ACE2

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Reactive

Human **Species**

Human

59272

Gene ID Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Antibody detection limit dilution 1:32000.

Limit

Western **Blot**

Approx 100kDa band was also observed on full-length recombinant Human ACE2 protein (calculated MW of 92.4kDa according to NP 001358344.1). Recommended concentration 1ug/ml. Primary incubation 1 hour at

room temperature.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa