



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ACHE ANTIBODY


SKU: EB06659

Telephone: (650) 697-3600



Email: customerservice@vectorlabs.com







Telephone: (650) 697-3600



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / acetylcholinesterase; apoptosis-related

Alias

Names

(YT blood group)|YT|ACHE

Usage

Immunofluorescence: Strong expression of the protein seen in the cytoplasm, plasma membrane and nucleus of U2OS cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of HeLa cells.

Summary

Recommended concentration: 10ug/ml.

Accession

NP_000656.1; NP_001269378.1

Blocking

EBP06659 **Peptide**

Peptide with sequence QFDHYSKQDRCSDL, from the C Terminus of the protein sequence according to NP 000656.1;

Immunogen NP 001269378.1.

Product

Comments

This antibody is expected to recognise isoform NP_000656.1 only (the ubiquitously expressed, hydrophillic form).

Peptide

Sequence

QFDHYSKQDRCSDL

Human, Rat, Mouse

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing

Method **Shipping**

Refrigerated Instructions

Predicted

Species

Reactive

Human Species

Human

43 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

ELISA

Detection Limit

Antibody detection limit dilution 1:32000.

Western

Approx 70kDa band observed in HepG2 and Jurkat cell lysates (calculated MW of 67.8kDa according to NP 000656.1).

Blot

Recommended concentration: 0.3-1µg/ml. Primary incubation 1 hour at room temperature.

Application

Type

Pep-ELISA, WB, IF, FC

GALLERY IMAGES





Telephone: (650) 697-3600



