

GOAT ANTI-ACK1 ANTIBODY

SKU: EB05180



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	ACK activated p21cdc42Hs kinase ACK1
Names	
Usage Summary	Has to be failed as SC (120kDa) and J. Biol. Chem, May 2001; 276: 18392 - 18398 (150kDa)
Accession ID	NP_005772
Immunogen	Peptide with sequence C-HLLGSWGPAHHKR, from the C Terminus of the protein sequence according to NP_005772.
Peptide Sequence	C-HLLGSWGPAHHKR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Mouse
Human Gene ID	10188
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 85-95kDa band seen in Mouse Heart extracts. Recommended for use at 4-6µg/ml. Overnight incubation is recommended with this antibody.