Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ACK1 ANTIBODY

SKU: EB05180

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias ACK|activated p21cdc42Hs kinase|ACK1

Usage

Names

Summary

Has to be failed as SC (120kDa) and J. Biol. Chem, May 2001; 276: 18392 - 18398 (150kDa)

Accession

NP 005772

Peptide with sequence C-HLLGSWGPAHHKR, from the C Terminus of the protein sequence according to **Immunogen**

NP 005772.

Peptide

C-HLLGSWGPAHHKR Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted Species

Human, Mouse

Reactive

Mouse **Species**

Human

10188

Gene ID **Product**

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

Western Western blot: Approx 85-95kDa band seen in Mouse Heart extracts. Recommended for use at 4-6µg/ml.

Blot Overnight incubation is recommended with this antibody.