



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ADAM15 ANTIBODY

SKU: EB05396

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Instructions

metalloproteinase-like, disintegrin-like, and cysteine-rich protein 15|metalloprotease RGD disintegrin protein|a Synonym / **Alias** disintegrin and metalloproteinase domain 15|metargidin|MDC15|a disintegrin and metalloproteinase domain

Names 15 (metargidin)|ADAM metallopeptidase domain 15 (metargidin)|ADAM15 NP_003806.3; NP_997078.1; NP_997079.1; NP_997080.1

Accession ID

Blocking

EBP05396

Peptide

Peptide with sequence C-APPPPTVSSLYL, from the C Terminus of the protein sequence according to **Immunogen**

NP_003806.3; NP_997078.1; NP_997079.1; NP_997080.1.

Product This antibody is expected to recognize isoforms 2, 4, 5, 6 (NP_003806.3; NP_997078.1; NP_997079.1;

Comments NP_997080.1 resp.)

Peptide

C-APPPPTVSSLYL Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Human **Species**

Reactive

Human

Species

Human 8751 Gene ID

Mouse

11490

Gene ID Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Western

Western blot: Approx 90-95kDa band seen in Daudi lysate. Recommended for use at 2-4µg/ml. **Blot**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600