

GOAT ANTI-ADAM15 ANTIBODY

SKU: EB05396



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	metalloproteinase-like, disintegrin-like, and cysteine-rich protein 15 metalloprotease RGD disintegrin protein a disintegrin and metalloproteinase domain 15 metargidin MDC15 a disintegrin and metalloproteinase domain 15 (metargidin) ADAM metalloproteinase domain 15 (metargidin) ADAM15
Accession ID	NP_003806.3; NP_997078.1; NP_997079.1; NP_997080.1
Blocking Peptide	EBP05396
Immunogen	Peptide with sequence C-APPPPTVSSLYL, from the C Terminus of the protein sequence according to NP_003806.3; NP_997078.1; NP_997079.1; NP_997080.1.
Product Comments	This antibody is expected to recognize isoforms 2, 4, 5, 6 (NP_003806.3; NP_997078.1; NP_997079.1; NP_997080.1 resp.)
Peptide Sequence	C-APPPPTVSSLYL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	8751
Mouse Gene ID	11490
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 90-95kDa band seen in Daudi lysate. Recommended for use at 2-4µg/ml.

