

GOAT ANTI-ADAMTSL2 ANTIBODY

SKU: EB09644



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	ADAMTSL2 ADAMTS-like 2 FLJ45164 KIAA0605 OTTHUMP00000064620 OTTHUMP00000162492
Accession ID	NP_055509.2
Blocking Peptide	EBP09644
Immunogen	Peptide with sequence CTHVTGNYRKGNH, from the internal region of the protein sequence according to NP_055509.2.
Peptide Sequence	CTHVTGNYRKGNH
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	9719
Mouse Gene ID	77794
Rat Gene ID	311827
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Preliminary experiments gave an approx 25kDa band in Human Kidney, Liver and Lung lysates after 0.1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 105kDa according to NP_055509.2. The 25kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)