

GOAT ANTI-ADIPONECTIN RECEPTOR 1 ANTIBODY

SKU: EB06718

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	TESBP1A FLJ42464 FLJ25385 CGI45 CGI-45 protein PAQR1 CGI-45 ACDCR1 HGNC:24040 ADIPOR1 adiponectin receptor 1
Accession ID	NP_057083.2; NP_001121159.1
Blocking Peptide	EBP06718
Immunogen	Peptide with sequence RYGLEGGCTDDTLL, from the C Terminus of the protein sequence according to NP_057083.2; NP_001121159.1.
Product Comments	Variants (NP_057083.2; NP_001121159.1) encode the same protein.
Peptide Sequence	RYGLEGGCTDDTLL
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Pig, Dog, Cow
Reactive Species	Human
Human Gene ID	51094
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 38kDa band observed in Human Skeletal Muscle lysates (calculated MW of 42.6kDa according to NP_057083.2). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

