

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ADIPONECTIN RECEPTOR 1 ANTIBODY

SKU: EB06718



Telephone: (650) 697-3600



250kDa 150kDa

100kDa

75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Unit Size

100 µg

Synonym / Alias

TESBP1A|FLJ42464|FLJ25385|CGI45|CGI-45 protein|PAQR1|CGI-45|ACDCR1|HGNC:24040|ADIPOR1|adiponectin

receptor 1

Names Accession

NP_057083.2; NP_001121159.1

Blocking

EBP06718 **Peptide**

Peptide with sequence RYGLEGGCTDDTLL, from the C Terminus of the protein sequence according to **Immunogen**

NP_057083.2; NP_001121159.1.

Product Comments

Variants (NP_057083.2; NP_001121159.1) encode the same protein.

Peptide

RYGLEGGCTDDTLL **Sequence**

Human

51094

Shipping

Instructions Refrigerated

Predicted

Species

Human, Mouse, Rat, Pig, Dog, Cow

Reactive

Species

Human

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Western Approx 38kDa band observed in Human Skeletal Muscle lysates (calculated MW of 42.6kDa according to

Blot NP_057083.2). Recommended concentration: 0.1-0.3μg/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa