

GOAT ANTI-ADRA1B (AA252-264) ANTIBODY

SKU: EB11920

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	alpha-1B-adrenergic receptor alpha-1B adrenoreceptor alpha-1B adrenoceptor alpha-1B adrenergic receptor OTTHUMP00000160844 ALPHA1BAR ADRA1 adrenergic, alpha-1B-, receptor ADRA1B
Accession ID	NP_000670.1
Blocking Peptide	EBP11920
Immunogen	Peptide with sequence C-TLRIHSKNFHEDT, from the internal region of the protein sequence according to NP_000670.1.
Peptide Sequence	C-TLRIHSKNFHEDT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Reactive Species	Human
Human Gene ID	147
Mouse Gene ID	11548
Rat Gene ID	24173
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 60kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 56.8kDa according to NP_000670.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

