

GOAT ANTI-AFAMIN / ALPHA ALBUMIN ANTIBODY

SKU: EB09538



250kDa
150kDa
100kDa
75kDa

50kDa
37kDa

25kDa
20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	alpha-albumin MGC125339 MGC125338 ALF ALBA ALB2 afamin AFM
Accession ID	NP_001124.1
Blocking Peptide	EBP09538
Immunogen	Peptide with sequence C-NKDDRPKDLSLRE, from the internal region of the protein sequence according to NP_001124.1.
Peptide Sequence	C-NKDDRPKDLSLRE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	173
Mouse Gene ID	75547
Rat Gene ID	293024
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 80kDa band observed in Human Brain (Cerebellum and Cerebral Cortex) lysates (calculated MW of 69.1kDa according to NP_001124.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

