


GOAT ANTI-AHI1 / JOUBERIN ANTIBODY

SKU: EB06539



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	contatins SH3 and WD40 domains OTTHUMP00000017265 OTTHUMP00000017264 OTTHUMP00000017263 Ahi-1 isoform I jouberin Abelson helper integration site dJ71N10.1 ORF1 JBTS3 FLJ20069 FLJ14023 DKFZp686J1653 AHI-1 AHI1
Usage Summary	Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.
Accession ID	NP_060121.3; NP_001128304.1; NP_001337433.1
Blocking Peptide	EBP06539
Immunogen	Peptide with sequence PTAESEAKVKTC, from the N Terminus of the protein sequence according to NP_060121.3; NP_001128304.1; NP_001337433.1.
Peptide Sequence	PTAESEAKVKTC

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	54806
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

