

GOAT ANTI-AICDA ANTIBODY

SKU: EB05005



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	activation induced cytidine deaminase HIGM2 CDA2 AID activation-induced cytidine deaminase AICDA
Names	
Accession ID	NP_065712
Immunogen	Peptide with sequence C-EVDDLRFDAFRTLGL, from the C Terminus of the protein sequence according to NP_065712.
Peptide Sequence	C-EVDDLRFDAFRTLGL
Shipping Instructions	Refrigerated
Predicted Species	Mouse
Reactive Species	Human
Human Gene ID	57379
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 25kDa band seen in Raji cell lysate. Recommended for use at 3-5µg/ml.