

GOAT ANTI-AKAP10 / D-AKAP2 ANTIBODY

SKU: EB06333



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	AKAP10 PRKA10 D-AKAP2 MGC:9414 A kinase (PRKA) anchor protein 10 dual specificity A-Kinase anchoring protein 2 A-kinase anchor protein 10 dual-specificity A-kinase anchoring protein 2 mitochondrial A kinase PPKA anchor protein 10 protein kinase A anchoring protein 10
Accession ID	NP_009133.2
Blocking Peptide	EBP06333
Immunogen	Peptide with sequence C-AQYDQPLEKSTKL, from the C Terminus of the protein sequence according to NP_009133.2.
Peptide Sequence	C-AQYDQPLEKSTKL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	11216
Mouse Gene ID	56697
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 65kDa band observed in Human Testis lysate (Calculated MW of 73.8kDa according to NP_009133.2) Recommended concentration: 0.3-1.0ug/ml
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

