

GOAT ANTI-AKAP3 / SOB1 ANTIBODY

SKU: EB06311



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	AKAP3 SOB1 FSP95 PRKA3 AKAP110 A kinase (PRKA) anchor protein 3 fibrousheathin 1 sperm oocyte binding protein 1 CT82 A-kinase anchor protein 3 A-kinase anchor protein, 110kDa Fibrous Sheath Protein of 95 kDa cancer/testis antigen 82 ibrous sheath protein, 95kDa protein kinase A anchoring protein 3
Accession ID	NP_006413.2
Blocking Peptide	EBP06311
Immunogen	Peptide with sequence C-TPLQLLDWLMVNL, from the C Terminus of the protein sequence according to NP_006413.2.
Peptide Sequence	C-TPLQLLDWLMVNL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Reactive Species	Human
Human Gene ID	10566
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 90kDa band observed in Human Testis lysates (calculated MW of 95kDa according to NP_006413). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

