Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ALP / PDLIM3 ANTIBODY

SKU: EB06186

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym / **Alias**

PDZ and LIM domain protein 3|actinin-associated LIM protein 3|enigma homolog|DKFZp686L0362|alphaactinin-2-associated LIM protein|ALP|PDZ and LIM domain 3|PDLIM3

Names

ID

Accession NP_055291.1; NP_001107579.1

Peptide with sequence C-PEGYDTVTLYPKA, from the C Terminus of the protein sequence according to **Immunogen**

NP 055291.1; NP 001107579.1.

Product

This antibody is expected to recognise both reported isoforms. Comments

Peptide

C-PEGYDTVTLYPKA Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Human, Mouse, Rat, Dog, Pig, Cow **Species**

Reactive

Human **Species**

Human

27295

Gene ID

Mouse 53318 Gene ID

Rat Gene ID 114108

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Western Blot Western Blot: Prominent band of approx 39 kD observed in Human Muscle lysate and a much less intense band of approx 35 kD in Human Heart lysate (predicted MWs of 42kDa according to NP_055291). Recommended for use at 0.1-0.5µg/ml.

Application Type

Pep-ELISA, WB