

GOAT ANTI-AMPA1 / GLUR1 ANTIBODY

SKU: EB07524



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	AMPA 1 MGC133252 HBGR1 GLURA GLUR1 GLUH1 glutamate receptor, ionotropic, AMPA 1 GRIA1
Accession ID	NP_000818.1
Blocking Peptide	EBP07524
Immunogen	Peptide with sequence C-EEGRDQTTSDQSN, from the internal region of the protein sequence according to NP_000818.1.
Peptide Sequence	C-EEGRDQTTSDQSN
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Human Gene ID	2890
Mouse Gene ID	14799
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments in Human Brain (Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

