


GOAT ANTI-ANGPT1 ANTIBODY

SKU: EB10271



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias ANGPT1|angiopoietin 1|ANG1|AGPT|AGP1

Names

Usage Summary **Immunofluorescence:** Strong expression of the protein seen in the plasma membranes of A431 cells. Recommended concentration: 10µg/ml. **Flow Cytometry:** Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.

Accession ID NP_001137.2; NP_001186788.1

Blocking Peptide EBP10271

Immunogen	Peptide with sequence C-TDQYNTNALQR, from the internal region of the protein sequence according to NP_001137.2; NP_001186788.1.
Peptide Sequence	C-TDQYNTNALQR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow, Pig
Reactive Species	Human
Human Gene ID	284
Mouse Gene ID	11600
Rat Gene ID	89807
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, IF, FC

GALLERY IMAGES

