
Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ANILLIN / SCRAPS (INTERNAL) ANTIBODY

SKU: EB07300



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names anillin (Drosophila Scraps homolog), actin binding protein|DKFZp779A055|anillin, actin binding protein (scraps

homolog, Drosophila)|ANLN|Scraps|ANILLIN

Usage Summary $<\!\!\textrm{strong}\!\!>\!\!\textrm{Immunofluorescence}\!\!:<\!\!/\textrm{strong}\!\!>\!\!\textrm{Strong expression of the protein seen in the nuclei of NIH3T3 and}$

Neuro-2a cells. Recommended concentration: 10µg/ml.

Accession

NP_061155.2

Blocking Peptide

EBP07300

Immunogen

Peptide with sequence C-TYPDDEKRKNP, from the internal region of the protein sequence according to

NP_061155.2.









Peptide C-TYPDDEKRKNP Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted Species

Human, Mouse, Rat, Dog

Reactive

Mouse **Species**

Human

54443

Gene ID Mouse

68743 Gene ID

Rat Gene ID 307056

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:16000.

Limit

Application

Pep-ELISA, IF Type

GALLERY IMAGES

