


GOAT ANTI-ANILLIN / SCRAPS (INTERNAL) ANTIBODY

SKU: EB07300



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names anillin (Drosophila Scraps homolog), actin binding protein|DKFZp779A055|anillin, actin binding protein (scraps homolog, Drosophila)|ANLN|Scraps|ANILLIN

Usage Summary **Immunofluorescence:** Strong expression of the protein seen in the nuclei of NIH3T3 and Neuro-2a cells. Recommended concentration: 10µg/ml.

Accession ID NP_061155.2

Blocking Peptide EBP07300

Immunogen Peptide with sequence C-TYPDDEKRKNP, from the internal region of the protein sequence according to NP_061155.2.

Peptide Sequence	C-TYPDDEKRKNP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Mouse
Human Gene ID	54443
Mouse Gene ID	68743
Rat Gene ID	307056
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, IF

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

