

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-ANNEXIN A11, BIOTINYLATED ANTIBODY**

**SKU:** EB05307-B



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

**Alias** 

ANXA11|annexin A11|ANX11|CAP50|56 kDa autoantigen|CAP-50|annexin XI|annexin-11|autoantigen, 56-

kD|calcyclin-associated annexin 50 Names

**Accession** 

NP\_001148.1

**Blocking** 

**Peptide** 

EBP05307-B

Peptide with sequence YPGYPPPGGYPPC., from the N Terminus of the protein sequence according to **Immunogen** 

NP 001148.1.

**Product** Reported variants represent identical protein: NP\_665875.1, NP\_665876.1, NP\_001265337.1, NP\_001265336.1,

NP\_001148.1 Comments

**Peptide** 

YPGYPPPPGGYPPC. Sequence

Human

311

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

**Species** 

Reactive

Human **Species** 

Human

Gene ID

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png Grade

**ELISA** 

**Detection** 

Antibody detection limit dilution 1:32000.

Limit

Blot

Western Approx 55kDa band observed in Jurkat lysates (calculated MW of 54.4kDa according to NP\_001148.1). See non-

biotinylated parental product's datasheet for further QC data. Recommended concentration: 1-3µg/ml.

**Application** 

**Type** 

Pep-ELISA, WB

## **GALLERY IMAGES**



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa