

GOAT ANTI-ANNEXIN A13 ANTIBODY

SKU: EB05304



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	MGC150460 annexin, intestine-specific annexin XIII ANX13 ISA annexin A13 ANXA13
Accession ID	NP_004297.2; NP_001003954.1
Blocking Peptide	EBP05304
Immunogen	Peptide with sequence GNRHAKASSPQGFD-C, from the N Terminus of the protein sequence according to NP_004297.2; NP_001003954.1.
Product Comments	This antibody is expected to recognise isoform a (NP_004297.2) and isoform b (NP_001003954.1).
Peptide Sequence	GNRHAKASSPQGFD-C
Shipping Instructions	Refrigerated
Predicted Species	Human, Pig
Human Gene ID	312
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Approx. 30kDa band observed in Human Kidney lysate at 0.5ug of primary. Please note that the band we observe is different from the predicted size of approx. 36kDa according to NP_004297. Currently we cannot find an explanation for this in the literature. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any splice variants/modified forms been reported?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)