

GOAT ANTI-ANTIGEN KI-67 / MKI67 (AA98-110) ANTIBODY

SKU: EB11653



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	proliferation-related Ki-67 antigen OTTHUMP00000020740 OTTHUMP00000020739 MKI67 KIA antigen KI-67 antigen identified by monoclonal antibody Ki-67
Accession ID	NP_002408.3; NP_001139438.1
Blocking Peptide	EBP11653
Immunogen	Peptide with sequence ENESLQNGRKSTE-C, from the N Terminus of the protein sequence according to NP_002408.3; NP_001139438.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_002408.3; NP_001139438.1).
Peptide Sequence	ENESLQNGRKSTE-C
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	4288
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Not yet tested - our routine western blotting protocol does not allow for the detection of proteins >250kDa. Therefore we cannot recommend an optimal concentration, and the antibody is an aspiring product.
Application Type	Pep-ELISA