


GOAT ANTI-AP3B1 / HPS ANTIBODY

SKU: EB06382



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	AP3B1 PE HPS HPS2 ADTB3 ADTB3A adaptor-related protein complex 3, beta 1 subunit beta3A-adaptin AP-3 complex beta-3A subunit
Usage Summary	Flow Cytometry: Flow cytometric analysis of HepG2 cells. Recommended concentration: 10ug/ml.
Accession ID	NP_003655
Blocking Peptide	EBP06382
Immunogen	Peptide with sequence SSNSFPYNEQSGG-C, from the N Terminus of the protein sequence according to NP_003655.
Peptide Sequence	SSNSFPYNEQSGG-C

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	8546
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

