


Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-AP3B1 / HPS ANTIBODY**





## **SPECIFICATIONS**

**Unit Size** 100 μg

Synonym / Alias

 $AP3B1|PE|HPS|HPS2|ADTB3|ADTB3A|adaptor-related\ protein\ complex\ 3,\ beta\ 1\ subunit|beta3A-adaptin|AP-3B-adaptin|AP-3B-adaptin|AP-3B-adaptin|AP-3B-adaptin|AP-3B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|AP-B-adaptin|$ 

complex beta-3A subunit

Names Usage

<strong>Flow Cytometry:</strong> Flow cytometric analysis of HepG2 cells. Recommended concentration:

Summary

10ug/ml.

Accession ID

NP\_003655

**Blocking** 

EBP06382

Peptide

**Immunogen** 

 $Peptide \ with \ sequence \ SSNSFPYNEQSGG-C, \ from \ the \ N \ Terminus \ of \ the \ protein \ sequence \ according \ to$ 

NP\_003655.

Peptide Sequence

SSNSFPYNEQSGG-C





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

Shipping Refrigerated Instructions

**Predicted** 

Human **Species** 

Reactive

Human

**Species** Human

8546

**Gene ID Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite\_medium.png$ 

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:64000.

Limit

Type

 $\begin{array}{c} \textbf{Application} \\ \textbf{Pep-ELISA, FC} \end{array}$ 

## **GALLERY IMAGES**

