


GOAT ANTI-APOA1 ANTIBODY

SKU: EB08746



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names preproapolipoprotein|apolipoprotein A1|apolipoprotein A-I, preproprotein|amyloidosis|MGC117399|apolipoprotein A-I|APOA1

Accession ID NP_000030.1

Blocking Peptide EBP08746

Immunogen Peptide with sequence C-SALEEYTKKLNTQ, from the C Terminus of the protein sequence according to NP_000030.1.

Peptide Sequence C-SALEEYTKKLNTQ

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human
Reactive Species	Human
Human Gene ID	335
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Liver shows sinusoid staining between the hepatocytes. Recommended concentration, 2-4µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

