


Telephone: (650) 697-3600

## **GOAT ANTI-APOA1 ANTIBODY**

**SKU:** EB08746



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

**Alias** 

ID

preproapolipoprotein|apolipoprotein A1|apolipoprotein A-I, preproprotein|amyloidosis|MGC117399|apolipoprotein A-I|APOA1

**Names** Accession

NP\_000030.1

**Blocking** 

EBP08746 **Peptide** 

**Immunogen** 

Peptide with sequence C-SALEEYTKKLNTQ, from the C Terminus of the protein sequence according to NP\_000030.1.

**Peptide** C-SALEEYTKKLNTQ

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

**Predicted** 

Human

**Species** Reactive

Human **Species** 

Human

335 Gene ID

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**IHC Results** 

In paraffin embedded Human Liver shows sinusoid staining between the hepatocytes. Recommended

concentration, 2-4µg/ml.

**ELISA** 

Detection

Antibody detection limit dilution 1:16000.

Limit

**Application** Pep-ELISA, IHC

Type

## **GALLERY IMAGES**

