



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-APOA5 ANTIBODY

SKU: EB07909



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

APOA5|apolipoprotein A-V|APOA-V|APOAV|MGC126836|MGC126838|RAP3|apolipoprotein A5|apolipoprotein

AV|regeneration-associated protein 3 Names

Accession

NP_443200.2 ID

Blocking

EBP07909 **Peptide**

Peptide with sequence C-HSLHDQGHSHLGDP, from the C Terminus of the protein sequence according to **Immunogen**

NP 443200.2.

Peptide Sequence

C-HSLHDQGHSHLGDP

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human **Species**

Human

116519 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Antibody detection limit dilution 1:64000. **Detection**

Limit

Western Approx 38kDa band observed in Human Liver lysates (calculated MW of 41.2kDa according to NP_443200.2).

Blot Recommended concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB Type

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa