

GOAT ANTI-APOA5 ANTIBODY

SKU: EB07909



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	APOA5 apolipoprotein A-V APOA-V APOAV MGC126836 MGC126838 RAP3 apolipoprotein A5 apolipoprotein AV regeneration-associated protein 3
Accession ID	NP_443200.2
Blocking Peptide	EBP07909
Immunogen	Peptide with sequence C-HSLHDQGHSHLGDP, from the C Terminus of the protein sequence according to NP_443200.2.
Peptide Sequence	C-HSLHDQGHSHLGDP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	116519
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 38kDa band observed in Human Liver lysates (calculated MW of 41.2kDa according to NP_443200.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

