


GOAT ANTI-APOBEC3D (AA61-74) ANTIBODY

SKU: EB10903



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	probable DNA dC->dU-editing enzyme APOBEC-3D apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3E pseudogene DNA dC->dU-editing enzyme APOBEC-3D ARP6 APOBEC3E APOBEC3DE apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3D APOBEC3D
Accession ID	NP_689639.2
Blocking Peptide	EBP10903
Immunogen	Peptide with sequence C-PKRQSNHRQEVYFR, from the internal region of the protein sequence according to NP_689639.2.
Peptide Sequence	C-PKRQSNHRQEVYFR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	140564
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Spleen. Recommended concentration: 3.75µg/ml.
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 45kDa band observed in lysates of cell line A549 (calculated MW of 46.6kDa according to NP_689639.2). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

