Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-APOBEC4 ANTIBODY

SKU: EB08273

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

putative C->U-editing enzyme APOBEC-4|OTTHUMP00000060680|MGC26594|C1orf169|apolipoprotein B mRNA **Alias**

editing enzyme, catalytic polypeptide-like 4 (putative)|APOBEC4 **Names**

Accession

NP 982279.1

Blocking

EBP08273 **Peptide**

Peptide with sequence C-QTKRNPNTKAQE, from the internal region of the protein sequence according to Immunogen

NP 982279.1.

Peptide

C-QTKRNPNTKAQE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

403314

Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western Preliminary testing showed a band at approx 40kDa in Human Testes lysate after 0.03µg/ml antibody staining Blot

(calculated MW of 41.6kDa according to NP_982279.1). Primary incubation 1 hour at room temperature.

Application

Pep-ELISA **Type**