

GOAT ANTI-APOH (AA145-157) ANTIBODY

SKU: EB11195



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names activated protein C-binding protein| anticardiolipin cofactor| APC inhibitor| apo-H| apolipoprotein H (beta-2-glycoprotein I)| B2G1| B2GP1| B2GPI| beta(2)GPI| beta-2-glycoprotein 1| BG|APOH

Accession ID NP_000033.2

Blocking Peptide EBP11195

Immunogen Peptide with sequence C-PSIPTFATLRVYK, from the internal region of the protein sequence according to NP_000033.2.

Peptide Sequence C-PSIPTFATLRVYK

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Reactive Species Human

Human Gene ID 350

Rat Gene ID 287774

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:1000.

Western Blot Approx 50kDa band observed in Human Liver, Kidney and Placenta lysates (calculated MW of 38.3kDa according to NP_000033.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Gu et al, Invest Ophthalmol Vis Sci. 1999 May;40(6):1229-37. PMID: 10235557). Recommended concentration: 0.2-0.6µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

