

GOAT ANTI-APOL2 ANTIBODY

SKU: EB08916



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	APOL2 apolipoprotein L, 2 APOL-II OTTHUMP00000028936 apolipoprotein L-II apolipoprotein L2
Accession ID	NP_112092.1; NP_663612.1
Blocking Peptide	EBP08916
Immunogen	Peptide with sequence C-TKIHEMLQPGQD, from the C Terminus of the protein sequence according to NP_112092.1; NP_663612.1.
Product Comments	Both reported isoforms represent identical proteins (NP_112092.1; NP_663612.1).
Peptide Sequence	C-TKIHEMLQPGQD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	23780
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 37.1kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 38kDa according to NP_112092.1; NP_663612.1). Recommended concentration: 0.03-0.1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

