


Telephone: (650) 697-3600

GOAT ANTI-APOLIPOPROTEIN C-III ANTIBODY

SKU: EB08795



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias

MGC150353|APOCIII|apolipoprotein C-III|APOC3

Names

Accession

NP_000031.1

Blocking

ID

Peptide

EBP08795

Peptide with sequence C-KDYWSTVKDKFSE, from the internal region of the protein sequence according to

Immunogen NP 000031.1.

Peptide Sequence

C-KDYWSTVKDKFSE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human

Species

Reactive **Species**

Human

Human

345 **Gene ID**

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results

In paraffin embedded Human Liver shows vesicular staining in the cytoplasm of hepatocytes. Recommended

concentration, 3-6μg/ml.

ELISA

Detection

Antibody detection limit dilution 1:4000.

Limit

Type

Application Pep-ELISA, IHC

GALLERY IMAGES

